

Engineering Consultancy
KODAK Motion Picture
Laboratory Process Survey

Catalog Number: 404 9771

Purpose

To provide the customer laboratory with an independent assessment of sensitometric conformance to Kodak recommended processing standards.

Recommended Procedures

Run the test in accordance with the documented instructions.

Benefits

- Compare your processing level to Kodak processing standards
- Receive recommendations for corrective actions required to bring processing Parameters back into control
- Cross checking of laboratory quality control measures

Kodak Deliverables

Kodak will provide a written report within 5 business days of receipt of processed survey strips, unless otherwise stated. The report includes full sensitometric curves, corrective recommendations, and verification of the laboratory in-house reference aim.

Additional services available:

KODAK Onsite Consulting Services, KODAK Photographic Solution Chemical Analysis can be ordered separately. Contact your local Kodak office for additional information.

Ordering Information

Contact your local Customer Service Representative or Kodak office for price and availability in your area.

For product and services updates, visit: www.kodak.com/go/motion/



Kodak Motion Picture Services

Kodak is a trademark. © Eastman Kodak Company 2004